


# gp140(HIV-1/N\_YBF30)



---

Product Name	gp140(HIV-1/N_YBF30)
Catalog No	IT-001-207p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	gp140(Group N Reference Strain: HIV-1 isolate YBF30) (aa 24-663; GenBank#O91086) with 6xHis tagged at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp140(HIV-1/N\_YBF30)

Note: