

# gp120(HIV-1/P\_RBF168)



---

Product Name	gp120(HIV-1/P_RBF168)
Catalog No	IT-001-208p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	gp120(Group P Reference Strain: HIV-1 isolate RBF168) (aa 38-515; GenBank#ACY40659) with a 6xHis tag at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

## Reference



E: gp120(HIV-1/P\_RBF168)

Note: