


small secreted glycoprotein [Bundibugyo virus]



Product Name	small secreted glycoprotein [Bundibugyo virus]
Catalog No	IT-014-016p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	small secreted glycoprotein [Bundibugyo virus] (aa 15-373, GenBank Accession No. QYU59239) with 6xHis tagged at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	≥95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

Reference



E: small secreted glycoprotein [Bundibugyo virus]

Note: